

Patient MRI Checklist

When a doctor requests an MRI scan:

- ✓ Tell the doctor who requested the MRI scan that you have an implanted spinal cord stimulation device.
- ✓ Make sure that you have your patient ID card and Remote Control. If you do not have your patient ID card, please call your Nevro HFX Care Team for assistance.
- ✓ **Before your MRI scan, all patients implanted with an HFX device must first undergo an impedance check.** Make sure you contact your Nevro HFX Care Team before scheduling your MRI scan. Your radiographer and Nevro HFX Care Team will confirm the results before your MRI.

When making your MRI appointment:

- ✓ Provide the information below to the MRI facility when making your appointment:
 - You have a Nevro HFX spinal cord stimulator
 - If you have any other medical device implants
 - A detailed description of your device including the product number, leads, and other implant components (this information can be found on your patient ID card)
 - The contact information of the doctor who implanted your device
 - The MRI guidelines, which can be found at www.Nevro.com/Manuals
- ✓ This information will be used by your physician and radiographer to determine if you can safely undergo an MRI scan.



PTRC1000




PTRC2300



PTRC2500



PTRC3000




Make sure your device is fully charged before your MRI scan.

Patient MRI Checklist

Preparing for your MRI appointment:

- ✓ Contact your Nevro HFX Care Team to confirm you can safely have the MRI scan scheduled.
- ✓ Bring your patient ID card and Remote Control to the MRI appointment.
- ✓ Before your MRI scan, be sure to charge your IPG and make sure your stimulation is switched to OFF. You can turn stimulation OFF with your Remote Control.
- ✓ After the MRI scan, turn your system back ON with your Remote Control and contact your doctor or Nevro HFX Care team to ensure that your device is programmed to pre-MRI settings.



Do you need help?

Please contact your local HFX Care Team who will help you find a center that is compatible with spinal cord stimulators, help gather the information you need, and who can answer your questions.



MRI Guidelines for Your Radiographer

Scan the code for detailed MRI compatibility information or visit [Nevro.com/manuals](https://nevro.com/manuals).

Note: HFX and HFX iQ are approved for an MRI scan under certain conditions. Different models of implanted devices, including lead types, result in different approved scanning conditions. Please review the Nevro HFX MRI Guidelines carefully with your physician, radiographer and Nevro HFX Care Team before getting an MRI scan.

SENZA, SENZA II, SENZA OMNIA, OMNIA, HF10, the HF10 logo, 10 kHz Therapy, HFX, the HFX logo, HFX iQ, HFX Algorithm, HFX CONNECT, the HFX Connect logo, HFX COACH, the HFX Coach logo, HFX CLOUD, the HFX Cloud logo, RELIEF MULTIPLIED, the X logo, NEVRO, and the NEVRO logo are trademarks or registered trademarks of Nevro Corp. ©2024 Nevro Corp. All Rights Reserved.