

HFX App

Compatibility List



Below are the Apple iPhones and supported OS version that are compatible with the HFX App:

- iPhone 8
- iPhone 8 Plus
- iPhone X
- iPhone XS
- iPhone XS Max
- iPhone XR
- iPhone 11
- iPhone 11 Pro
- iPhone 11 Pro Max
- iPhone SE 2nd Gen
- iPhone 12 Mini
- iPhone 12
- iPhone 12 Pro
- iPhone 12 Pro Max
- iPhone 13 Mini*
- iPhone 13
- iPhone 13 Pro
- iPhone 13 Pro Max
- iPhone SE 3rd Gen
- iPhone 14
- iPhone 14 Pro
- iPhone 14 Plus
- iPhone 14 Pro Max
- iPhone 15
- iPhone 15 Pro
- iPhone 15 Plus
- iPhone 15 Pro Max
- iPhone 16
- iPhone 16 Pro
- iPhone 16 Plus
- iPhone 16 Pro Max
- iPhone 16e
- iPhone 17
- iPhone 17 Air
- iPhone 17 Pro
- iPhone 17 Pro Max



Compatible iOS Versions

16, 17, 18, 26

Minimum Available Memory (MB)

- 500



Additional Information

- If your device isn't listed on this page, it's not supported by the HFX App.
- The HFX App is not compatible with iPad or iPod devices.
- * If you're using an iPhone 13 Mini, please keep your phone within 2 meters (6 ½ feet) from your HFX iQ device when using the HFX App. This will make sure the connection between the two devices is strong enough so your system works properly.



Nevro Corp.
1800 Bridge Parkway Redwood City, CA 94065
1.650.251.0005 • 1.888.956.3876

SENZA, SENZA II, SENZA OMNIA, OMNIA, HF10, the HF10 logo, 10 kHz Therapy, HFX, the HFX logo, HFX iQ, the HFX iQ logo, HFX Algorithm, HFX CONNECT, the HFX Connect logo, HFX ACCESS, the HFX Access logo, HFX COACH, the HFX Coach logo, HFX CLOUD, the HFX Cloud logo, RELIEF MULTIPLIED, the X logo, NEVRO, and the NEVRO logo are trademarks or registered trademarks of Nevro Corp.

©2026 Nevro Corp. All Rights Reserved.

ML-00031 Rev.I S000899 V4.1



NevroHFX.com